

HORIZON[®]

Delivering versatility, performance and reliability

HORIZON[®] is a versatile, user friendly TMS system with a comprehensive range of services designed to help treat more patients safely and effectively each day, as well as run a successful TMS practice. HORIZON[®] delivers standard and theta burst protocols and the energy recovery system prevents pulse decay during stimulation to provide reliability from pulse 1 to pulse 3000.



HORIZON[®] Air Film Coil

Capable of delivering repetitive pulses at high power levels, the coil features a thermostatic control for uninterrupted protocols.



HORIZON[®] MT Remote Coil

An ergonomic hand-held coil for motor threshold assessment with power controls located directly on the coil.

Coil Stand

The multi-movement articulated arm provides precise bi-lateral coil positioning. Attaching the coil stand to the chair creates a smaller foot-print system.

Head Support System

Comfortable and supportive headrest with contra-lateral side pad to ensure coil-to-patient contact throughout treatment.

Therapy Chair

Mobile and motorized, the chair is fully adjustable for patient comfort



HORIZON[®] Stimulator

High performance stimulator delivering 60Hz standard and 50Hz theta burst protocols for shorter treatment sessions. The stimulator energy recovery system presents pulse decay during stimulating sessions.

Stimulator Cart

Mobile cart with locking swivel wheels. Holds the stimulator at the ideal height for stimulation.




TMS Patient Cap

Ideal for treatment location marking, allowing for repeated coil positioning, shortening the treatment set-up time.



HORIZON[®] AFC Cooling Module

Powers the fan of the Air Film Coil when not connected to the stimulator, reducing the time between treatments.



ORDERING INFORMATION

HORIZON® FAST PACKAGE - P/N: 5195-00	
Description	Part Number
HORIZON® Mainframe	4950-00
HORIZON® Power Supply Unit (PSU)	4806-00
HORIZON® User Interface (UI)	5048-00
HORIZON® Air Film Coil (AFC)	4980-00
HORIZON® MT Remote Coil	4990-00
HORIZON® Cooling Module	5018-00
Foot Switch	9525-01
Therapy Chair (220V / 110V)	280105/280106
Chair Coil Stand	280111
Stimulator Cart	3767-00
US Mains Filter / EU Mains Filter	5158-00/5176-00

PATIENT CAPS	
Description	Part Number
Patient cap - size S*	4941-10
Patient cap - size M*	4942-10
Patient cap - size L*	4943-10
Patient cap - size XL*	4944-10
EASYCAP (EEG Landmark cap) - Size S	315900
EASYCAP (EEG Landmark cap) - Size M	315901
EASYCAP (EEG Landmark cap) - Size L	315902

* Please note: patient caps have a minimum order quantity of 10 caps.

CHAIR ACCESSORIES	
Description	Part Number
Chair leg rest plastic cover	280110
Chair paper cover (roll)	280108
Chair paper roll holder	280109



FOR MORE INFORMATION; VISIT WWW.MAGSTIM.COM/CONTACT, OR CALL US ON: 01994 240798

The Magstim Company Ltd.
Spring Gardens, Whitland, Carmarthenshire, SA34 0HR
MK1720-02



HORIZON® is a CE marked medical product, indicated for the stimulation of the human cortex for therapeutic and research purposes, including the treatment of depression.



www.magstim.com

Not for US audiences

magstim
THE BRAINS BEHIND TMS™